

## 3B BIRTHING SIMULATOR P90 PRO

Obstetric simulation with the **3B Scientific Birthing Simulator P90 PRO** to enhance training of manual maneuvers throughout the birthing process. Includes a fully articulating fetus with palpable fontanelles and open mouth.

- Effective training of manual birthing maneuvers
- Demonstration of all possible fetus presentations and positions
- Realistic anatomy including landmarks such as spine and pelvic bones
- Includes a fully articulating fetus
- Inflatable uterus for intrauterine manipulations and Leopold's maneuver



### 3B Scientific Birthing Simulator P90 PRO

#### features:

- Effective training of manual birthing maneuvers
- Demonstration of all possible fetus presentations and positions
- Realistic anatomy including landmarks such as spine and pelvic bones
- Includes a fully articulating fetus
- Cesarean incision insert
- Inflatable uterus for intrauterine manipulations and Leopold's maneuver
- Inserts for training of artificial rupture of membranes (AROM/ amniotomy)

Using the 3B Scientific Birthing Simulator P90 PRO, trainees can see how manual maneuvers affect the fetus in-utero and train the following manual delivery skills and techniques:


- Normal vaginal delivery
- Normal delivery of umbilical cord and placenta
- Palpation and assessment of fetal position and presentation
- Demonstration of all normal and abnormal fetal positions and presentations (like face, neck, or brow) to represent routine and difficult deliveries
- Intrauterine manipulation of fetus
- Shoulder dystocia (adaptive birth canal and pelvic bones) demonstration and management
- Breech deliveries (including complete, footling or incomplete, and frank breech presentation)
- Palpation of fetal fontanelles features
- Complete and incomplete placenta birth
- Demonstration of placenta previa: total, partial, and marginal
- Prolapse of umbilical cord
- Nuchal umbilical cord
- Cesarean delivery
- Artificial rupture of membranes (ARM or AROM, amniotomy)
- Clamping and cutting of umbilical cord
- Demonstrate suctioning of baby's mouth and nose
- Manual maneuvers (such as Leopold, Pinard's, Mauriceau, Prague, Woods and Rubin maneuver)

### 3B Scientific Birthing Simulator P90 PRO 1022879

- Detachable, padded stomach cover with advanced insert for cesarean delivery
- Inflatable uterus insert for Leopold maneuver including cover that can be cut for cesarean delivery
- Birth canal (portio, cervix, vulva)
- Fastening belt
- Lubricant
- Amniotic fluid and 100 inserts
- Assembly kit
- Fully articulating fetal baby
- Umbilical cord and placenta

### 3B Scientific Birthing Simulator P90 BASIC 1022878

With pre-cut cesarean insert and basic birth canal, but without amniotomy insert or inflatable uterus.



*Realistic pelvic anatomy to demonstrate shoulder dystocia and practice manual maneuvers*

03/2022



ALL CONSUMABLES AND UPGRADE KITS ARE AVAILABLE AT [3BSCIENTIFIC.COM](http://3BSCIENTIFIC.COM)



3B Scientific GmbH • Ludwig-Erhard-Straße 20 • 20459 Hamburg • Germany • [3bscientific.com](http://3bscientific.com)  
Phone: +49 40 73966-0 • Fax: +49 40 73966-100 • E-mail: [info@3bscientific.com](mailto:info@3bscientific.com)