



Be  
**READY**

DANi Teen & Young Adult Advanced Training Solution

*Nasco*  
HEALTHCARE

# Quick Start Guide

DANi™ is packaged with arms and female genitalia attached. Legs have been disconnected to prevent possible damage during shipping.

1. Carefully unpack the product's individual components (see Fig. 1). Your DANi manikin includes:

- Head, torso, arms and female genitalia
- Left leg
- Right leg
- Male genitalia
- Accessories



Figure 1



2. Connect left and right legs to the torso (see Fig. 2).

3. Female genitalia attached, or connect male genitalia (see Fig 3 & 4).

4. Your product is ready for use.

DANi comes pre-assembled to support nasogastric tube procedures (lavage, gavage). Please see section Oral and Nasal Lavage of the instruction manual if you want to configure your product to support different airway management procedures.

For more information on your DANi manikin, please contact our customer service team at 1-833-NASCOHC (627-2642) or email us at [info@nascohealthcare.com](mailto:info@nascohealthcare.com). We are happy to help you!



Figure 2



Figure 3



Figure 4



Please visit  
<https://nhcpartnerportal.com/products/manuals/>  
for complete instructions.



Nasco Healthcare  
16 Simulaid Drive  
Saugerties, NY 12477  
1-833-NASCOHC (627-2642)  
[info@nascohealthcare.com](mailto:info@nascohealthcare.com)  
[www.nascohealthcare.com](http://www.nascohealthcare.com)

©2022 Nasco Healthcare

PP9201

PC 101037 02/23